CK

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

- The analyzer will finish the self check after start-up.
- It is recommend to do analyzer calibration monthly and for each new lot of test kit.

Step 1 Sample Preparation

1a Allow the test kit back to room temperature for 30 minutes before use.

1b Use the R2 with quantitative capillary to collect sample.


1c Insert the R2 into R1.

Note:

- Please update the standard curve with the barcode on the R1 cuvette if a new lot test kit is to be used.
- The capillary of the R2 should be fully filled.

Step 2 Testing

2a Insert the R1 into test channel of analyzer.

2b The analyzer will mix the sample from R2 capillary with R1 automatically.

2c The analyzer will mix the R2 and R1 automatically.

2d The analyzer will test and print the results automatically.

Calibration

The calibration values for the different lots of the kits are stored on the two-dimensional code on the cuvette. Before test the new lot of kits, the instrument automatically scan the two-dimensional code on the cup to obtain the corresponding calibration curve during testing.

Quality control

3-level calibration system guarantee the results' reliability for each lot of test kits, including the instrument calibration, remote reagent calibration and the third party calibration.

The third party calibration applicable for:

1. The daily indoor quality control test.
2. New lots of reagent.

3. New operator training.
4. The results can not match the clinical symptoms.
5. The first use of the reagent.

If still can not be calibrated, contact the manufacture for further technical support.

Reference Value

<164U/L

Recommended that each laboratory establish its own reference range.

Interpretation

1. If the test result is $\geq 164\text{U/L}$, indicating that there may be a risk of myocardial infarction or muscle damage, and it is recommended that the patient undergo further examination.
2. The result only for clinical reference, comprehensive consideration should be combined with the clinical management of patients with symptoms / signs, medical history, other laboratory tests and treatment response.
3. All laboratory tests depend on random errors. If the test results are in doubt, or if they do not match the clinical symptoms, re-test the sample or confirm the results with other methods.

Limitations

Bilirubin $\leq 340\mu\text{mol/L}$, Hemoglobin $\leq 5\text{ g/L}$, Triglyceride $\leq 5\text{mmol/L}$ will not affect the test result.

Performance Characteristics

1. Linearity range: $4\text{U/L} \sim 1200\text{U/L}$, $r \geq 0.990$.
2. Within-run precision: relative standard deviation (RSD) $\leq 5\%$.
3. Between-run precision: Inter-batch relative range (R) $\leq 10\%$.
4. Accuracy: relative deviation (B) should not exceed $\pm 10\%$.
5. Blank limit: reagent blank absorbance ≤ 0.5 .
6. Reagent blank absorbance change rate ≤ 0.002 .
7. Analytical sensitivity: when the sample concentration is 500U/L , its absorbance change rate is between 0.060 and 0.160.

Precautions

Attention:

Only for in vitro diagnostic.

Only for professional use.

All samples and reactive wastes are treated as sources of infection.


Do not use the kits beyond shelf life.

Do not mix different batches of reagents.

Warning:

To avoid error, do not forced to take out the cuvette from the device. Follow the device operation manual strictly, If the problem cannot be solved, contact the manufacturer for further technical support.

SYMBOLS USED ON LABELS

Symbol	Usage	Symbol	Usage
	Use-By date		Do not freeze
	Batch code		Biological risks
	Manufacturer		Do Not Reuse
	Temperature Limit		
	Contains sufficient for <n> tests		
	Do not use if package is damaged		
	Consult Instructions for use		
	Keep Away from Sunlight		
	In Vitro Diagnostic Medical device		
	Authorized Representative in the European Community		

References

Shang Hong, Wang Liu three, Shen Ziyu and so on. National Clinical Laboratory Practice (Fourth Edition) [M] Beijing: People's Medical Publishing House, 2015:412-413.

Approval Date&Revision Date

Approval Date: Feb.15,2022

Revision Date: Jan 1, 2023